



stryker®

Medical

Your Partner **In Delivery**

LD304

The labor and delivery bed
that delivers rewarding
birth experiences.

The Stryker LD304 Birthing Bed

– Comfort and safety for mom, baby and you!

Having a baby is one of the most important moments in a woman's life. As a nurse, one of your greatest concerns is ensuring that mom has a rewarding birth experience.

That's why Stryker designed the LD304. Tailor-made for moms-to-be, the Stryker birthing bed not only offers the highest level of patient safety and comfort, but also gives caregivers an added sense of security through enhanced ergonomics, simple operation and advanced mobility.

With the Stryker bed as your partner, you can stay focused on what matters most – helping moms bring healthy little ones into the world.

Product Highlights:

1. Intuitive

The Stryker birthing bed's innovative design and functionality helps ensure the safety of you and your patients during every phase of labor, delivery and transport.

2. Ergonomic & Safe

Nearly infinite leg positioning capabilities on the LD304 birthing bed allow you to comfortably accommodate moms of all sizes. And, at the same time, these positioning features allow you to stay close to your patients – and keep your hands free – so you can focus on caring for mom throughout her delivery.

3. Comfortable

Our premium mattress system provides five inches of plush foam, along with a built-in egg crate surface. The minimal gap between the head and foot mattresses also lends to a feeling of comfort and support for mom.

Intuitive

"Stryker has put together a bed that's really a no-brainer in terms of ease of use."

In one easy motion, leg supports are in position without interference from the mattress. To customize to mom, simply slide the supports toward or away from her using the easy-squeeze handle. If mom prefers, place her leg in the attached calf supports and adjust to the desired angle – lock into place with a flip of the lever.





Remove the foot section with ease... and save your back.

Removing the foot section is simple and safe. With its ergonomic design that promotes proper body mechanics, the foot section takes less effort to slide off. You'll appreciate that this helps protect you from potential back strain and injury during its removal and replacement. When the foot section is on the bed, you can feel assured that it's securely attached, as Stryker's exclusive Lock-Rite™ feature provides a visual indicator that it is locked into place.



Ergonomic & Safe

"The foot section is so easy to take off...and better yet, so much easier to put back on. And you know it's locked when the bar is down."

Stryker's removable foot section solves many of the challenges L & D nurses face every day – it's easier to position moms, it's more responsive to the varying needs of the delivery process and leaves no doubt that it's properly secured to the bed.

Comfortable

"The more comfortable you can make a patient in labor, the longer she can go without having to use alternative methods, such as medicines."

The LD304's inflatable, full-width lumbar pillow gives mom proper support for back labor, and the ample seat pillow enhances perineal visibility. You can also engage Stryker's powered Skoocher™ feature to move mom closer to the perineal edge, and provide her with full and proper back support.

Suite additions for lasting impressions

Did you know that the industry leader in patient handling equipment also offers you a complete line of furniture? Built with the same high standards you've come to expect from Stryker, these suite additions create lasting impressions of quality in any labor and delivery room.



Upright Sitting

Footrests or calf supports can be positioned for the patient, with the foot section still intact. Allows nurse to provide hands-free care from the foot end of the bed.



Squatting

By supporting herself on the labor bar mom can fully stretch the pelvis and vaginal opening, allowing the baby to descend.



Resting/Squatting

Use of high semi-fowler and Skoocher allows for comfortable resting position.



McRobert's Maneuver

For expected shoulder dystocia, the one-touch adjustable seat and length-adjustable footrests facilitate pelvic opening. Nurse's hands are free to apply suprapubic pressure.



Lithotomy

Calf supports and Skoocher help position mom for better visibility and access – whether for an episiotomy, forcep application or if mom is anesthetized.



Side Lying

Calf support is readily accessible and easy to position for appropriate flexion. Holds mom's leg to increase pelvic stretch, while helping to reduce the nurse's back strain.



Knee/Chest

Footrest can be lowered to increase the back stretch, to help relieve back labor and aid in fetal rotation.



Kneeling, Leaning Forward

Engaging Skoocher and upright backrest puts mom in increased pelvic tilt, allowing for better fetal descent and rotation in posterior position.

Created to support the best

in maternity care.



Intuitive controls, including one-touch Trendelenburg button.

In intermediate height, Stryker's three-position siderails help mom enter and exit the bed, placing her correctly and comfortably in the seat section. This mid-range siderail height also helps mom feel secure in the bed, and still feel close to her baby and family.



Joint Replacements

Trauma

Spine

Micro Implants

Orthobiologics

Instruments

Interventional Pain

Navigation

Endoscopy

Communications

Patient Handling Equipment

EMS Equipment

6300 South Sprinkle Road
Kalamazoo, Mich. 49001 USA
t: 269 329 2100 f: 269 329 2311
toll free: 800 787 9537

www.stryker.com

The information presented in this brochure is intended to demonstrate the breadth of Stryker product offerings. Always refer to the package insert, product label and/or user instructions before using any Stryker product. Products may not be available in all markets. Product availability is subject to the regulatory or medical practices that govern individual markets. Please contact your Stryker representative if you have questions about the availability of Stryker products in your area.

Products referenced with ™ designation are trademarks of Stryker.
Products referenced with ® designation are registered trademarks of Stryker.

Literature Number: LLB-456-01A
PC/SCG-25M-0504

Copyright © 2004 Stryker
Printed in USA