

Gynnie®

OB/GYN Stretcher



Adjustable footrests feature simple, one-handed operation

4" thick Ultra Comfort mattress provides unparalleled support

Sealed mattress and clean lines prevent fluid collection

Compact design with pop-up push handles and retractable fifth wheel improves maneuverability

Intuitive Patient Positioning

Gynnie combines the trademark durability and features of a standard Stryker stretcher with supreme OB/GYN functionality. Ideal for pelvic exams and general use, this multi-functional stretcher can outfit an entire ER department.



Glide footrests out, position patient and lower foot section to provide excellent pelvic access.



Gynnie can be configured with optional OB/GYN accessories.



Medical

3800 E. Centre Ave.
Portage, MI 49002
t: 269 329 2100 f: 269 329 2311
toll free: 800 787 9537

www.stryker.com

Standard Features

- 4" (10 cm) Ultra Comfort Mattress with sealed seams
- 8" (20 cm) Omni Surface™ casters with wheel covers
- 26" (66 cm) wide patient surface
- Dual cylinder pneumatic backrest
- Dual pedestal hydraulics with constant rate descent
- Dual side-mounted foot control hydraulics
- Retractable fifth-wheel steering
- Four-wheel, steel-ring brakes with dual side and head activators
- Glideaway™ siderails
- Integrated storage tray with oxygen bottle holder
- Permanently-attached folding foot section
- Permanently-attached foot supports
- Pop-up push handles
- Roller bumpers
- Three IV receptacles
- Unique one-handed foot section release lever

Specifications

• Model Number	1061	
• Overall Length	81.5"	(207 cm)
• Overall Width		
Siderails Up	31.5"	(80 cm)
Siderails Down	26.75"	(68 cm)
• Weight Capacity	500 lbs	(228 kg)
• Height Range (to litter top)		
High	36"	(91 cm)
Low	22.5"	(57 cm)
• Litter Positioning		
Backrest	0° - 90°	
Foot Section	0° - 90°	
Trend./Reverse Trend.	±18°	
• Patient Surface	26" x 77"	(66 x 196 cm)
• Siderails	13" x 39"	(34 x 99 cm)
• Caster Diameter	8"	(20 cm)

Stryker reserves the right to change specifications without notice.

Warranty

- One year parts, labor and travel or two years parts

Optional Features

- C-spine cassette film holder
- Calf supports and brackets
- Backrest X-ray cassette holder
- Hoop for drainage basin
- IV poles
 - Permanent folding two- or three-stage
 - Standard removable
- Patient restraints
- Serving/instrument tray
- Siderail pads
- Upright oxygen bottle holder
- Extended warranty available