

Vortex Mixers

Mini Vortex Mixer



◆ Description

The mini vortex mixer is the ideal vortex mixer for every lab with its compact size, ease of use, and its economic price point. The 3,000 rpm fixed speed ensures there is no sacrificing any mixing power.

◆ Features

- Extremely compact and easy operation;
- Powerful vortexing for tubes up to 30mm;
- Head is made of durable material to ensure long life;
- Refined steel base is designed to ensure product stability;
- Choice of different colored working head.

◆ Specifications

Product Name	Mini Vortex Mixer			
Model	VM-03U	VM-03RU	VM-03GU	VM-03PU
Color	Blue	Red	Green	Purple
Oscillation Mode	Circumference Oscillation			
Orbital Diameter	Ø4.5mm			
Capacity	Max.80ml			
Tube Diameter	≤Ø30			
Ambient Temperature	5~40°C			
Speed	3000rpm			
Power Supply	AC100~240V/50~60Hz			
Maximum Power	5W			
Outside Dimensions(ΦxH)	96X85mm			
Weight	0.55kg			
IP rating (DIN EN 60529)	IP40			