

Mini Centrifuges



- ◆ Description
Our mini centrifuges are essential laboratory equipment used in clinical, environmental and industrial research. These personal centrifuges are ideal for quick spins, microfiltration of samples, cell separations and many other routine laboratory procedures.
- ◆ Features
 - Ergonomically designed with a small footprint;
 - High grade plastic cover for clear visibility;
 - Smooth and quiet operation;
 - Automatic stop when lid is opened;
 - Digital control delivers precise speed;
 - Fits a variety of tubes with included accessories.

◆ Accessories

Name	Part Number	Picture	Pcs/Pack
Standard Rotor for 6 x 1.5/2.0ml Tubes	MLXF-7		1
Rotor for 16 x 0.2ml Tubes	MLXF-3		1
Adaptors for 0.5ml Tubes	MLXF-1		6
Adaptors for 0.2ml Tubes	MLXF-2		6



MLX-104



MLX-106



MLX-110

◆ Specifications

Model	MLX-104	MLX-106	MLX-110
Speed	4000rpm	6000rpm	10000rpm
G-force(RCF)	900xg	2000xg	5000xg
Noise Level	≤65dB	≤65dB	≤75dB
Body Color	Pink	Off White	Blue
Ambient Temperature	5°C~40°C		
Operational Humidity	≤80%		
Power Supply	AC110V/60Hz(US) AC220V/50Hz(Eur)		
Weight	≤1Kg		
Standard Heads	Head for 8 Tubes (in 2 rows): 16x0.2ml Standard Head: 6x1.5ml		
Dimensions WxDxH	155x140x115mm		
Safety Feature	Stops When Lid is Opened		