

MH6616

Slide Drying Bench

The Slide Drying Bench aids in the preparation of microscope slides at the specimen mounting stage. It complements both the Paraffin Wax Dispenser and the Paraffin Wax Mounting Bath to provide a complete set of paraffin wax sampling equipment for the modern laboratory.

Accepting up to 50 slides (76 x 25mm), the Slide Drying Bench has the facility for drying slides in different orientations- across the drying support bars, angled from the bars, or flat without the bars.

There is an On/Off switch, and "Power On" white neon and "Heater On" amber neon indicators.

The silicone rubber mat element heater provides heating up to 100°C, and the temperature is controlled from 10°C to 100°C by a built-in energy regulator.

The case and top are powder-coated aluminium and the unit comes complete with handy carrying handles.

Technical Information

Material	Powder coated aluminium base and top
Max Element Temp.	100°C
Dimensions (d x w x h), mm	180 x 390 x 95
Weight, kg	1.8

Key Features

- Up to 50 slides capacity
- Different drying orientations
- Operates up to 100°C max
- Replacement elements are easily fitted



MH6616

Ordering Information

Model	Capacity	Electrical Requirement
MH6616	50 slides	230V 50/60Hz, 150W
MH6616X1	50 slides	115V 50/60Hz, 150W
MH6616X6	50 slides	230V 50/60Hz, 150W, EU Plug