

Stretcher Surface Portfolio



	Pioneer® Pressure Redistribution Mattress	Ultra Comfort Mattress	Enhanced Comfort Mattress
Description	The Pioneer® mattress offers the greatest pressure redistribution in the Stryker stretcher mattress portfolio. The eleven horizontal air bladders use self-adjusting technology to redistribute the patient's weight, helping to relieve pressure.	The Ultra Comfort mattress combines dual density foam with a blue, Dartex® polyurethane cover.	Standard on many Stryker stretchers, the Enhanced Comfort mattress has dual density foam with a Vinyl cover.
Composition	Low ILD 10 foam at the head section, 11 horizontal air bladders, and 15° foot end slope.	Low ILD 15 foam	Dual density foam
Cover	Dartex® Cover with 4-way stretch to help prevent shearing	Dartex® Cover with 4-way stretch to help prevent shearing	Vinyl Cover
Weight Capacity	700 lb. weight capacity	700 lb. weight capacity	700 lb. weight capacity
Thickness	5" Thickness	4" or 5" Thickness	2", 3", or 4" Thickness
Width	26" or 30" Width	26" or 30" Width	26" or 30" Width