

VTs 098

-86°C ULTRA LOW TEMP. FREEZER

VTs 098 combines extreme performance and exceptional use of space with a fantastic and functional design. The -86°C unit features a one-compressor system, internal LED-lights and 100% stainless steel compartments. The characteristic Vestfrost door-handle keeps the door perfectly tight and the top-panel lock enables users to avoid unauthorized access.



SPECIFICATIONS

ISO 9001, 14001 & OHSAS 18001 CERTIFIED

Gross volume, litres (cu. ft.)	93.3 (3.3)
Net volume, litres (cu. ft.)	91 (3.2)
Temperature range, at 16°C to 22°C AMB	-60°C to -86°C
Refrigerant	Nature R
Internal stainless steel shelves	2

DIMENSIONS

Height, mm (inches)	830 (32.7)
Width, mm (inches)	600 (23.6)
Depth, mm (inches)	643 (25.3)
Depth incl. door handle, mm (inches)	700 (27.6)
Gross weight, kg (lbs)	65 (143.3)
Net weight, kg (lbs)	56 (123.5)

LOADING QUANTITIES

Qty. per 20' / 40' container	54 / 108
------------------------------	----------

Products and specifications are subject to change due to continuous product development.



BIOMEDICAL SOLUTIONS

VESTFROST
SOLUTIONS

A/S Vestfrost
Falkevej 12
DK-6705 Esbjerg Ø
www.vestfrostsolutions.com